

# **INCREASING THE EFFICIENCY OF FUSIONS FOR XRF**


**NAVAS INSTRUMENTS – JOSE LAS NAVAS**



# USING FUSION FOR XRF ANALYSIS IN THE CEMENT INDUSTRY


Due to complex chemistry of Portland and Pozzolanic Cements dissolution of the clinker phases necessary to get meaningful XRF results

ONLY WAY TO ELIMINATE OR COMPENSATE FOR THIS EFFECT IS FUSION




Different x-ray intensities for the same concentrations of Ca in different Ca bearing compounds in cement


# CEMENT WITH LITHIUM BROMIDE



# NIST SAMPLE



# MOISTURE AT 600 °C





# THE NEED FOR FUSION

## Advantages:


- Very homogeneous sample
- No mineral phases or particles
- Inter-element matrix effects reduced by dilution
- Any remaining inter- element matrix effects easily corrected
- Easy to handle and store

**Essentially the ideal XRF sample**



# CHARACTERISTICS OF A GOOD FUSED BEAD

- All solid material completely dissolved into glass
- No loss of sample, flux or additives during fusion – **EXACT** ratio is **KNOWN**
- No residue left behind in crucible
- Comes out of the casting mold or moldable crucible cleanly



# FUSION MAKERS-DESIGN CONSIDERATION

## Precise Temperature Control

Advantages of circular VERSUS linear designs for multiple sample prep stations in how you measure the temperature

## Sufficient agitation for good mixing

Rotation as well as horizontal and vertical movements – ideal figure of eight movement

## Atmosphere control – sufficient air for complete oxidation

Inject air to keep uniform temperature and air circulation via convection



# UNIQUE ASPECTS OF THE NAVAS FUSION SYSTEM

- Electric furnace, NO danger of explosion, No gas
- Analyzes loss on ignition and provides beads in the same instrument, simultaneously
- Simple, Reliable, MODULAR, Robust construction
- Bench type instrument with high throughput at a reasonable price
- No pouring, only moldables used
- Automation, requires little operator involvement
- Minimal moldable cleaning needed, only after 20 - 30 uses, not required after each use
- Easy to use software in Windows
- Up to 8 samples simultaneously



*Navas Instruments AFS-5000-6,  
Automatic Fusion & LOI System For  
XRF Spectrometry*

# FACTORS IN FUSION

## Weighing

*Weighing error*



Sample



Flux (Li<sub>2</sub>B<sub>4</sub>O<sub>7</sub>)



*GOI*

## Fusion

*Flux evaporation*

FeO  
↓  
Fe<sub>2</sub>O<sub>3</sub>

O<sub>2</sub>




CO<sub>2</sub>


H<sub>2</sub>O

*LOI*

Pt crucible




Fused bead




1000-1150 °C

# CONTROLLING FLUX, SAMPLE, OXIDANT AND NON-WETTING AGENT RATIOS



Dilution ratio: weight ratio of flux to sample

# TEMPERATURE VS LOSS OF FLUX, ADDITIVES AND COMPONENTS



# HOW TO CONTROL LOSSES

- **Use Pre-fused fluxes**
- **Monitor weight in-situ in furnace**
- **Add releasing agent just before fusion is completed**

**Reduces loss of material associated with LiBr evaporation**

- **Measure Fused Bead Weight**
- **Use Fused Bead Weight to Sample Weight to account for LOI, GOI and LOF(Loss of Fusion)**

# DILUTION RATIO CORRECTION

Calibration equation with dilution ratio correction

$$W_i = (b I_i + c) \left( 1 + \sum_{j \neq L} \alpha_j W_j + \alpha_F R_F + K_F \right)$$

$$K_F = -\alpha_F \bar{R}_F$$

General calibration equation

$$W_i = (b I_i + c) \left( 1 + \sum \alpha_j W_j + \alpha_F \Delta R_F \right)$$

$\Delta R_F = R_F - \bar{R}_F$

$\left. \begin{array}{l} \Delta R_F \\ R_F \end{array} \right\} \begin{array}{l} \text{Difference between the} \\ \text{actual and standard} \\ \text{ratio} \\ \text{Actual dilution ratio} \end{array}$

$\alpha_F R_F + K_F$  is the correction term for the difference between the actual and standard dilution ratio.

# HOW USING THIS APPROACH IMPROVES YOUR ACCURACY

10 replicate beads read on XRF

	Na2O	MgO	Al2O3	SiO2	P2O5	SO3	K2O	CaO	TiO2	MnO	Fe2O3
	mass%	mass%	mass%	mass%	mass%	mass%	mass%	mass%	mass%	mass%	mass%
Using actual sample weight and bead weight											
Std dev.	0.024	0.023	0.016	0.051	0.009	0.008	0.007	0.086	0.007	0.002	0.012
<b>RSD(%)</b>	<b>13.9</b>	<b>1.48</b>	<b>0.28</b>	<b>0.25</b>	<b>4.44</b>	<b>0.25</b>	<b>1.25</b>	<b>0.13</b>	<b>2.41</b>	<b>2.63</b>	<b>0.28</b>
Using sample weight and flux weight											
Std,Dev	0.036	0.023	0.018	0.079	0.010	0.011	0.008	0.086	0.008	0.004	0.013
<b>RSD(%)</b>	<b>15.92</b>	<b>1.58</b>	<b>0.33</b>	<b>0.42</b>	<b>5.23</b>	<b>0.37</b>	<b>1.46</b>	<b>0.15</b>	<b>2.75</b>	<b>5.01</b>	<b>0.31</b>

# SUMMARY

**The Navas Automated Fusion System has distinct advantages over current automated fusion machines available in the market place**

- The ability to record and update weights of the sample, flux and fused bead in situ allows for accurate tracking of the fusion process. Not available in other instruments
- Possible to do accurate LOI measurements at the same time as fusing the sample into a bead
- The use of fused bead weights in the final XRF calculation greatly improves the accuracy and precision of data used to control critical processes in cement manufacture



# CONCLUSIONS

**This novel approach that allows LOI and Fusion in the same process**

- Saves time and thereby improves productivity which translates into financial savings
- The ability to improve both the accuracy and long term precision of the measurements, improves the product reliability which also translates into cost savings and intangible assets with respect to product reliability

**THANK YOU FOR  
YOUR ATTENTION**